

Public Awareness Campaign to Combat Coronavirus (COVID-19)



Following the emergence of the Omicron variant Updates have been made to the COVID-19 Rapid Test Protocol to include swabs for both the throat and nose




Follow the instructions provided with the COVID-19 rapid testing device

How to use the device?




Step 1

Swipe the soft tip of the swab across the very back of the throat behind the Arch by opening mouth wide and stick out the tongue, ensure to swab the sides of the arch. Don't swab the tongue or teeth.




Step 2


Using the same swab, place tip of the cotton swab inside nostril and gently wipe around inside of the nose, no force needed & repeat it to other nostril.



The result is negative if one line appears in area (C).



The result is positive if two lines appear in the (C) / (T) areas. If so, a PCR test must performed.



Note: The test shall be considered **invalid** and must be repeated if no line appears in area (C).

Note




You should not be eating or drinking 30 mins before the procedure to yield reliable result.




The used tools must be sealed in a plastic bag and kept outside for 3 days. After 3 days, it can be treated as normal waste.







Instructions on testing positive following a rapid COVID-19 test





If you have tested positive

following a rapid COVID-19 test



You must report the result via the "BeAware" application





- •Muharraq Drive-Through Testing Centre, near the King Hamad University Hospital
- Southern Governorate Drive-Through Testing Centre, near the Bahrain Centre for Strategic, International, and Energy Studies (Derasat)
- Northern Governorate Drive-Through Testing Centre in Salman City



Without the need to call 444



Do not forget to carry the rapid test device with you following placement in a transparent airtight bag

Note:

Please isolate yourself, avoid interactions with others and wear a mask.